

## Specimen/Assay Description

Question Text	Answers Text
Collection by attentive prostate massage:	Yes, No, Not provided
Specimen:	Urine – sedimented, Urine – unsedimented, Prostatic ejaculate
Name of transport media:	
Transport time (hours):	
Holding temperature (Celsius):	
Storage temperature (-C):	
Method used:	Quantitative Real Time PCR, Nucleic acid sequence based amplification, Transcription-mediated Amplification
Assay used identified as:	1 <sup>st</sup> , 2 <sup>nd</sup> or 3 <sup>rd</sup> generation test uPM3, PCA3Plus, Progenesa, Aptima, “Gen-Probe test”, Test not specified
Housekeeping gene:	PSA mRNA, Other, Not provided
Result unit:	PCA3 Score, Other
Testing blinded to outcomes:	Yes, No, Not provided
Number and percent of samples with insufficient mRNA:	
Number and percent of reportable results:	
PCA3 losses explained?	Yes, No
ABSTRACTOR COMMENTS: If you would like to leave a comment pertaining to the information above indicate your name below:	